

Upgrade to a **TA.XT^{plus}** Texture Analyser

Upgrade your TA.XT2i Texture Analyser to a state of the art **TA.XT^{plus}** and take full advantage of the new performance specifications offered by the improved electronics and 32-bit **Texture Exponent** software.

At the same time, protect the initial capital investment of your original texture analyser purchase.



New Features


Develop your own test routines

A standard facility of the **TA.XT^{plus}** is the ability for the user to develop and write their own unique test procedure or "Sequence". The number of these sequences that can be produced and their complexity is limitless. What ever you want the texture analyser to do is now possible without being limited by standard library tests. Designing your own test method is made extremely easy and user friendly by using the Sequence Wizard provided in the 32-bit software.

Improved accuracy and speed

The new technical specifications of the **TA.XT^{plus}** will allow more accurate analysis of products at all sorts of speeds. Test speeds can be varied between 0.01mm/s up to 20mm/s with greater distance resolution. Data recorded during product testing can now be run at 500pps for longer timespans allowing for more precise analysis.

TA.XT^{plus} Texture Analyser



Reduced maintenance costs

New electronics allow accurate reporting of errors and the ability to update the firmware via flash upgrading from the internet. The software is also still constantly being improved and free updates will be available for downloading from our website. In addition upgrading will extend your instrument's warranty by an additional year.

32-bit Texture Exponent Software

Exponent is a completely new 32-bit software package exclusively designed for use with the **TA.XT^{plus}**. Some of the extra features available in addition to everything currently available in Exceed include:

- Extended File Names
- Enhanced Calibration Procedures
- Improved Presentation, Layout and Report Options
- New Improved Application Help Files
- Compatibility with existing Graph files
- One cable connection to a Laptop or Desktop PC

For information on the **TA.XT^{plus}** and **Texture Exponent** and their many new features please contact your local sales team. For details please consult our website: www.stablemicrosystems.com or contact us at the address below.

Stable Micro Systems

Vienna Court, Lammas Road, Godalming, Surrey GU7 1YL, UK

Telephone: *44 1483 427345

Facsimile: *44 1483 427600

e-mail: sales@stablemicrosystems.com

web: <http://www.stablemicrosystems.com>

