

UPGRADE TO A

TA.XT^{plus}
Texture Analyser®

Upgrade your TA.XT2i Texture Analyser to a state-of-the-art **TA.XT^{plus}** and take full advantage of the new performance specifications offered by the improved electronics and 32-bit **Exponent** software.

At the same time, protect the initial capital investment of your original texture analyser purchase.

New Features

Develop your own test routines

A standard facility of the TA.XT^{plus} is the ability for the user to develop and write their own unique test procedure or 'Sequence'. The number of these sequences that can be produced and their complexity is limitless. Whatever you want the texture analyser to do is now possible without being limited by standard library tests. Designing your own test method is made extremely easy and user friendly by using the Sequence Wizard provided in the 32-bit software.

Improved accuracy and speed

The new technical specifications of the TA.XT^{plus} will allow more accurate analysis of products at all sorts of speeds. Test speeds can be varied between 0.01mm/s up to 20mm/s with greater distance resolution. Data recorded during product testing can now be run at 500pps for longer timespans allowing for more precise analysis.





Reduced maintenance costs




New electronics allow accurate reporting of errors and the ability to update the firmware via flash upgrading from the internet. The software is also still constantly being improved and free updates are available for downloading from our website. In addition, upgrading will extend your instrument's warranty by an additional year.

32-bit Exponent software

Exponent is a 32-bit software package exclusively designed for use with the TA.XT^{plus}.

Some of the extra features available in addition to everything currently available in *Exceed* include:

- Extended file names
- Enhanced Calibration procedures
- Improved Presentation, Layout and Report options
- New improved Application Help files
- Compatibility with existing Graph files
- One cable connection to a laptop or desktop PC



EXPONENT
Stable Micro Systems

EXPONENT
Lite

Stable Micro Systems constantly improve the capability of *Exponent* software by adding testing procedures, presentation formats and analytical solutions, this being of particular importance for those in Research. However, many users require simple test and data analysis solutions for QA or in a high throughput production environment.

Therefore, two versions are now available – *Exponent* and *Exponent Lite* – both of which offer the same basic functions, although *Exponent* has additional data collection and analysis features and enhanced security options.

For information on the TA.XT^{plus} and *Exponent*, and their many new features, please contact your local sales team. For details please consult our website www.stablemicrosystems.com or contact us at the address below.

Stable Micro Systems®

Vienna Court, Lammas Road, Godalming, Surrey GU7 1YL, UK

Tel: *44 1483 427345

Fax: *44 1483 427600

Email: sales@stablemicrosystems.com

app.support@stablemicrosystems.com

tech.support@stablemicrosystems.com

Web: <http://www.stablemicrosystems.com>

